

Fruit acids of lemon

Sample: lemon juice

Inj. solvent: mobile phase

Inj. volume: 10 µL

SP: 150x4.6 Reprosil-Pur C18-AQ (3 µm)

MP: 0.5% aq. H₃PO₄

Flow: 2 mL/min

Detection:


— 205 nm

— 260 nm

Detector cell volume: 10 µL

Temperature: 25 °C

Pressure: 190 bar



ELSICO HPLC Company



SCAN

Consulting and
training for HPLC

www.hplc.ru/training



<https://sites.google.com/site/scansychev/>
<http://www.chromforum.ru/>