



Hydroxycinnamic acids in different specimens

Sample:

- cofe, water extract
- echinacea, ethanol extract
- red wine

Inj. solvent: mobile phase

Inj. volume: 10 uL

SP: 50x4.6 Reprosil-Pur Basic C18 (1.9 um)

MP: 17% AcN in 0.5% aq. H3PO4

Flow: 2 mL/min

Detection: 330 nm

Detector cell volume: 10 uL

Temperature: 25 °C

Pressure: 190 bar

ELSICO HPLC Company



SCAN

Consulting and training for HPLC

www.hplc.ru/training



<https://sites.google.com/site/scansychev/>
<http://www.chromforum.ru/>