



Catechins of green tee

Sample: green tea (SenCha), ethanol extract
Inj. solvent: mobile phase
Inj. volume: 10 uL
SP: 50x4.6 Reprisil-Pur Basic C18 (1.9 um)
MP: 17% AcN in 0.5% aq. H3PO4
Flow: 2 mL/min
Detection:
 — 205 nm
 — 270 nm
Detector cell volume: 10 uL
Temperature: 25 °C
Pressure: 200 bar

ELSICO HPLC Company



SCAN

Consulting and training for HPLC

www.hplc.ru/training



<https://sites.google.com/site/scansychev/>
<http://www.chromforum.ru/>