



ONE OF THE WORLD'S LEADING HPLC COLUMNS

Develop or Improve your HPLC Method

- Wide pH stability from 1.5 to 10*
- Exceptionally rugged USP Phases: Silica, C5, C8, C18, Phenyl-Hexyl, CN, NH₂, SCX
- Easy method scalability with 3, 5, 10, 10-PREP and 15μ medias
- Extensive batch traceability and reproducibility data supplied with every column

If Luna does not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, send in your comparative data within 45 days and keep the Luna column for FREE!

Material Characteristics

Packing Material	Particle Shape/Size (μ)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Load %	Calculated Bonded Phase Coverage (umole/m ²)
Luna Silica(2)	Spher. 3, 5,10,15	100	400	0	—
Luna C5	Spher. 5,10	100	440	12.5	7.85
Luna C8	Spher. 5,10	100	440	14.75	5.50
Luna C8(2)	Spher. 3, 5,10,15	100	400	13.5	5.50
Luna C18	Spher. 5,10	100	440	19	3.00
Luna C18(2)	Spher. 3, 5,10,15	100	400	17.5	3.00
Luna Phenyl-Hexyl	Spher. 3, 5,10, 15	100	400	17.5	4.00
Luna CN*	Spher. 3, 5, 10	100	400	7.0	3.80
Luna NH ₂ *	Spher. 3, 5, 10	100	400	9.5	5.80
Luna SCX*	Spher. 5,10	100	400	0.15 meq/g	

USP Phases for Virtually Every Application

Luna Bonded Phase Selectivity Chart

	Description	Application
C5	5, 10μ A reversed-phase chemistry that offers greater hydrolytic stability and alternative selectivity over the classic C4 phase.	Highly hydrophobic compounds. Good alternative to C8 or C18 when less retention is desired.
C8(2)	3, 5, 10, 10-PREP, 15μ C8 phase for excellent efficiency, peak shape and resolution. Significantly improved performance over traditional C8 phases due to high surface coverage. Slightly lower carbon load than original Luna C8.	Great starting point for methods development of pharmaceuticals, nucleotides and polar compounds. Excellent for highly aqueous and LC/MS applications.
C8	5, 10μ Original Luna C8 phase. Highly retentive C8 optimized for long lifetimes at extremely high pH levels.	A "workhorse" for extremely high pH applications, or when a more retentive C8 is desired.
C18(2)	3, 5, 10, 10-PREP,15μ C18 phase. Excellent efficiency, peak shape and resolution. Slightly lower carbon load than original Luna C18.	A C18 phase for virtually all HPLC applications. Range of particle sizes offers amazing versatility for capillary and LC/MS, to prep and process scale applications. Our No.1 recommended choice for HPLC method development.
C18	5, 10μ Original Luna C18 phase. Very high retention, optimized for long lifetimes at extremely high pH levels.	A "workhorse" for extremely high pH applications, or when longer retention is desired.
Phenyl-Hexyl	3, 5, 10, 10-PREP,15μ A (patented) phenyl phase which employs a hexyl alkyl linker as opposed to the traditional propyl chain. Offers great stability as well as alternative selectivity.	Remarkable selectivity for aromatic compounds. Offers selectivity advantages of both a phenyl phase and a C6 or C8 alkyl phase. Stable, reproducible alternative to traditional phenyl phases. Also very good for amine and polar compound selectivity.
CN	3, 5, 10μ Cyano phase. Can be used as reversed- or normal-phase material. The use of Luna base silica results in overall phase reproducibility and performance.	Excellent for improving the retention of polar compounds. Extremely rapid equilibration makes it very good for rapid screening and gradient applications. Exceptional normal phase performance.
NH ₂	3, 5, 10μ Amino Phase. Can be used in reversed or normal phase modes. Stable from pH 1.5 to 11.0 and under 100% aqueous conditions. High performance silica and bonding techniques produces a rugged, highly reproducible column.	Excellent for sugar analysis or anionic compounds by reversed phase chromatography and for compounds capable of hydrogen bonding under normal phase conditions.
SCX	5, 10μ A Benzene Sulfonic Acid bonded phase is used to make this Strong Cation Exchange (SCX) column. Offers great peak shape and resolution.	Excellent for cationic and nitrogen containing analytes such as cough and cold compounds, codeine and morphine alkaloids, and drugs of abuse. Also useful for 2-dimension chromatography (SCX-RP-MS-MS) of small peptides and protein digests.

*Luna CN is stable from pH 1.5 to 7.0, Luna Amino (NH₂) is stable from pH 1.5 to 11.0, Luna SCX is stable from pH 2.0 to 7.0.



LUNA SILICA

A Backbone & Phase Designed for Long Column Lifetimes

Luna's excellent silica is not simply the result of ultra-pure metal-free silica (99.999% purity). Meticulous care is given to the quality control of surface smoothness, pore structure and pore consistency to ensure particles of uniform structure and enhanced mechanical strength. Either bonded or unbonded, Luna silica produces highly advanced HPLC columns:


- Low percentage of "fines" from damaged silica particles stabilizing the column bed and leading to lower backpressures and enhanced column performance and lifetimes
- High column bed stability enhanced by particle shape uniformity

Incredible Silica Smoothness


Luna silica is extremely smooth and spherical. For bonded phases, this provides a uniform bonding surface for consistent and even bonded phase coverage. The likelihood of silica particle shearing and breakage during bonding and packing is very low; thus, Luna columns have high efficiencies and long column lifetimes.

	Low	High
Hydrophobicity	N/A	
H-Bonding	N/A	
Aromatic Selectivity	N/A	
Silanol Activity (pH 2.5)		
see p. 365 for testing probes.		
pH Stability	—	
Particle Size	3 μ , 5 μ , 10 μ , 10 μ -PREP, and 15 μ	
Phase	Unbonded silica	
Application	Non-polar compounds	


Compare Surface Smoothness



Phenomenex Luna®
5μ C18



Waters Symmetry®
5μ C18




Agilent Technologies Zorbax®
5μ SB-C18

App ID 1138

Natural Products

Column: Luna 5 μ Silica(2)
Dimensions: 150 x 4.6mm
Order No: 00F-4274-E0
Mobile Phase: Hexane/Dioxane (85:15)
Flow Rate: 1.5mL/min
Detection: UV @ 230nm
Sample:
1. Dihydrokawain
2. Yangonin
3. Kawain
4. Dihydromethysticin
5. Methysticin




LUNA C18(2), C8(2), C5

Take One Small Step and a Giant Leap in Performance

Luna has found a place as one of the World's top reversed phase columns because it provides a measurable improvement over many HPLC columns for two important chromatographic properties: resolution and peak shape. The high efficiencies and bonded phase surface coverage provide for sharp peaks. The result:




- Free exposed silanols virtually eliminated by complete bonding and endcapping
- Sharp peak shape for good method sensitivity
- pH stable from 1.5 to 10.0 for over 10000 hours

Peak Asymmetry Comparison of Competing Columns**




Conditions for all columns

App ID 1100
 Dimensions: 150 x 4.6mm
 Mobile Phase: Acetonitrile/Water (50:50)
 Flow rate: 1.0mL/min
 Detection: UV @ 254nm
 Sample: 1. Pyridine
 2. Phenol
 3. Toluene



Pyridine Peak Asymmetry Comparison**



Comparison of 7 different 5 μ reversed phase columns. This survey measures the degree of silanol activity on the surface of each silica. In this survey, Luna 5 μ C18(2) material demonstrates the lowest silanol activity.

** The comparative data presented here may not be representative for all applications.




Luna C18(2), C8(2), C5 (continued)


Wide pH Stability for Long Column Lifetime and Method Flexibility

Luna was one of the first columns to provide an extended pH range of 1.5 to 10** for most of its phases. The broad pH stability range is partly achieved by the density of the bonded phase. Peak shapes remain sharp independent of pH, indicating free silanols are well shielded from the analyte compounds. With Luna it is possible to work above the pK_b of many basic compounds, with the compound in a “neutral charge state”, thus the basic compound becomes more hydrophobic (less polar). This will improve retention, peak shape, and reproducibility.


pH Stability Data for Luna 5 μ C18(2)




Test Conditions: Column stability tested under highly acidic conditions. Continuous flush in 0.1% TFA (pH 1.5) in Water/Acetonitrile, 50:50.



Test Conditions: Column stability tested under highly basic conditions. Continuous flush in 20mM Na₂HPO₄ (pH 10.0) in Water/Acetonitrile, 50:50.



Optimizing resolution and selectivity with mobile phase pH increases method development flexibility



Conditions for Separations

Column: Luna 5 μ C18

Dimensions: 150 x 4.6mm

Order No.: 00F-4041-E0

Mobile Phase: Buffer/ Acetonitrile, (95:5) pH as shown

Flow Rate: 1.0mL/min

Temperature: 40°C

Detection: UV @ 270nm

Sample: 1. Sulfanilic Acid 3. Sulfisoxazole

2. Sulfamerazine 4. Sulfapyridine



** Luna CN is stable from pH 1.5 to 7.0 - Luna NH₂ is stable from pH 1.5 to 11.0 - Luna SCX is stable from pH 2.0 to 7.0.




* UV-Vis spectral evidence indicates that Sulfisoxazole is unstable at low pH.

Luna C18(2), C8(2), C5 (continued)

Steroids

App ID 1185
Column: Luna 5 μ C18(2)
Dimensions: 150 x 4.6mm
Order No.: 00F-4252-E0
Mobile Phase: 0.1% H₃PO₄/Acetonitrile/Methanol (54:35:11)
Flow Rate: 0.75 mL/min
Detection: UV @ 254nm
Sample:


1. Hydrocortisone	3. 11- α -Hydroxyprogesterone
2. Corticosterone	4. Cortisone Acetate
	5. 11-Ketoprogesterone



Antibacterials

App ID 1201
Column: Luna 5 μ Phenyl-Hexyl
Dimensions: 150 x 4.6mm
Order No.: 00F-4257-E0
Mobile Phase: A: 20mM KH₂PO₄, pH 2.5 B: Acetonitrile
Gradient: A/B (80:20) to A/B (75:25) in 5 min, to A/B (55:45) in 15 min
Flow Rate: 1.0 mL/min
Detection: UV @ 254nm
Sample:


1. Carbadox	5. Sulfadimethoxine
2. Thiamphenicol	6. Sulfaquinoxaline
3. Furazolidone	7. Nalidixic Acid
	8. Piromidic Acid



Water Soluble Vitamins


App ID 3305
Column: Luna 5 μ C18(2)
Dimensions: 150 x 4.6mm
Order No.: 00F-4252-E0
Mobile Phase: 20mM Potassium Phosphate, pH 3.0/Acetonitrile (95:5)
Flow Rate: 1.5mL/min
Temperature: 22°C
Detection: UV @ 214 nm
Sample:

1. Thiamine	4. Pantothenic acid
2. Cyanocobalamin (Vitamin B12)	5. Niacinamide
3. Ascorbic acid	6. <i>p</i> -Aminobenzoic acid



Acetaminophen, USP Method


App ID 1293
Column: Luna 5 μ C18(2)
Dimensions: 150 x 4.6mm
Order No.: 00F-4252-E0
Mobile Phase: Water/Methanol/Acetic Acid (69:28:3)
Flow Rate: 1.5mL/min
Detection: UV @ 275nm
Sample:



Pharmaceutical Preservatives


App ID 1141
Column: Luna 5 μ C5
Dimensions: 150 x 4.6mm
Order No.: 00F-4043-E0
Mobile Phase: A:0.5% Acetic acid in water/acetonitrile (80:20)
B: 0.5% Acetic acid in water/acetonitrile (20:80)
Gradient: A/B (100:0) to A/B (0:100) in 30min
Flow Rate: 1mL/min
Temperature: 25°C
Detection: UV @ 254 nm
Sample:

1. Propylparaben impurity	6. 8
2. Benzyl alcohol	
3. Phenol	
4. Benzoic acid	
5. Methylparaben	
6. Benzaldehyde	
7. Ethylparaben	
8. Propylparaben	



Penicillin Antibiotics

App ID 1215
Column: Luna 5 μ C8
Dimensions: 150 x 4.6mm
Order No.: 00F-4040-E0
Mobile Phase: Methanol/25mM KH₂PO₄, pH 4.6 (55:45)
Flow Rate: 1mL/min
Temperature: 22°C
Detection: UV @ 254 nm
Sample:






Luna C18(2), C8(2), C5 (continued)

App ID 1239 | Sun Screens

Column: Luna 5 μ C18
Dimensions: 250 x 4.6mm
Order No: 00G-4041-E0
Mobile Phase: Methanol/Water (80:20)
Flow Rate: 1.5mL/min
Temperature: 22°C
Detection: UV @ 240nm
Sample:


- 1. Benzophenone 3
- 2. DCHP
- 3. Parsol MCX
- 4. Homosalate 2
- 5. Octyl salicylate
- 6. Homosalate 1
- 7. Homosalate



App ID 9903 | Hormones

Column: Luna 5 μ C18
Dimensions: 100 x 4.6mm
Order No: 00D-4041-E0
Mobile Phase: Water/Tetrahydrofuran/Acetonitrile (60:22.5:17.5)
Flow Rate: 1.2mL/min
Temperature: 40°C
Detection: UV @ 244nm
Sample:


1. Norgestrel	5. syn-Norgestimate
2. syn-Norgestrel oxime	6. syn-Norgestrel oxime acetate
3. anti-Norgestrel oxime	7. anti-Norgestimate
4. Norgestrel acetate	8. anti-Norgestrel oxime diacetate



App ID 1288 | Polyamines (Dansylated)

Column: Luna 5 μ C18 + Security Guard
Dimensions: 250 x 4.6mm
Order No: 00G-4252-E0 and KJ0-4282 and AJ0-4287
Mobile Phase: A: 93% Methanol in Water
B: 65% Acetonitrile in Water
Gradient: A/B (0:100) to A/B (100:0) in 28 min to A/B (0:100) in 40 minutes
Temperature: 35°C
Flow Rate: 1.5mL/min
Detection: 360nm Excitation, 530nm Emission
Sample:


1. Putrescine	5. Spermine
2. DHX	6. Homospermine
3. DHP	7. Methylmethyhomospermine
4. Spermidine	8. Diethylhomospermine (DEHOP)



App ID 1284 | Glucocorticoids

Column: Luna 5 μ C8
Dimensions: 150 x 4.6mm
Order No: 00F-4040-E0
Mobile Phase: Acetonitrile/Water (50:50)
Flow Rate: 1.0mL/min
Detection: UV @ 240nm
Sample:


- 1. Dexamethasone 21-Phosphate
- 2. Prednisolone
- 3. Dexamethasone
- 4. Prednisolone 21-Acetate
- 5. Betamethasone 17-Valerate



App ID 2578 | Cough and Cold Compounds

Column: Luna 5 μ C8(2)
Dimensions: 150 x 4.6mm
Order No.: 00F-4249-E0
Mobile Phase: A: Methanol B: Acetonitrile w/0.1% H₃PO₄
C: Water w/0.1% H₃PO₄ and 0.1% heptane sulfonate
A/B/C (40:10:50)
Flow Rate: 1.5mL/min
Temperature: 22°C
Detection: UV@ 214nm
Sample:


1. Acetaminophen	3. Benzoic acid
2. Pseudoephedrine	4. Chlorpheniramine
	5. Dextromethorphan



App ID 1072 | Saw Palmetto Berry, p-Bromophenyl Esters

Column: Luna 3 μ C8(2)
Dimensions: 150 x 4.6mm
Order No.: 00F-4248-E0
Mobile Phase: Acetonitrile/Water (87:13)
Flow Rate: 1.5mL/min
Temperature: 25°C
Detection: UV @ 254nm
Sample:

1. Capric Acid	5. Linoleic Acid
2. Lauric Acid	6. Palmitic Acid
3. Linolenic Acid	7. Oleic Acid
4. Myristic Acid	8. Stearic Acid




LUNA PHENYL-HEXYL

Engineered For Stability

Luna Phenyl-Hexyl can provide separations not achievable on C18 or C8 columns; such as increased retention for polar, aromatic compounds as well as reversals in analyte elution order. Luna Phenyl-Hexyl is a reproducible, extremely stable phenyl phase. Most phenyl phases use a short propyl (3 carbon) linker, which limits phase stability. The Phenyl-Hexyl bonded phase employs a phenyl ring with a hexyl (6 carbon) linker and is densely bonded to Luna silica surface. Dense bonding and the hexyl linker reduce bonded phase hydrolysis and increases chemical stability. The result:





- Highly reproducible and stable phenyl phase
- Dual selectivity of both phenyl phase and a short alkyl phase (C5 or C8)
- Excellent retention of aromatic and polar, amine compounds
- 1.5 to 10 pH stability for over 6000 hours



HPLC

Luna

Chromatographic Comparisons of Phenyl Columns*



Antibacterials

Conditions for all columns:

Dimensions: 150 x 4.6mm

Mobile Phase: A: 20mM KH₂PO₄, pH 2.5

B: Acetonitrile

Gradient: Start A/B (80:20) to A/B (75:25) in 5 min, go to A/B (55:45) in 15 min

Flow Rate: 1.0mL/min

Detection: UV @ 254nm




Sample:

- | | |
|------------------|---------------------|
| 1. Carbadox | 5. Sulfadimethoxine |
| 2. Thiamphenicol | 6. Sulfaquinoxaline |
| 3. Furazolidone | 7. Nalidixic Acid |
| 4. Oxolinic Acid | 8. Piromidic Acid |

*The comparative data presented here may not be representative for all applications.

Luna Phenyl-Hexyl (continued)


pH Stability of Luna Phenyl-Hexyl



App ID 1199

Food Additives


Column:	Luna 5 μ Phenyl-Hexyl	
Dimensions:	150 x 4.6mm	
Order No:	00F-4257-E0	
Mobile Phase:	A: 50mM KH ₂ PO ₄ + 0.1%H ₃ PO ₄ B: Acetonitrile	
Gradient:	A/B (75:25) to A/B (25:75) in 18min, hold at A/B (25:75) for 12min	
Flow Rate:	1mL/min	
Temperature:	22°C	
Injection:	20 μ L	
Detection:	UV @ 230nm	
Sample:	1. Saccharin 2. p-Hydroxybenzoic acid 3. Sorbic acid 4. p-Hydroxybenzoic acid methyl ester 5. Dehydroacetic acid 6. p-Tolnic acid 7. p-Hydroxybenzoic acid ethyl ester 8. n-Propyl p-hydroxybenzoate	



App ID 2572

Cough and Cold Medicine

Column:	Luna 3 μ Phenyl-Hexyl	
Dimensions:	75 x 4.6mm	
Order No:	00C-4256-E0	
Mobile Phase:	A: Acetonitrile B: Methanol/20 mM KH ₂ PO ₄ (80:20) pH 9.0	
Gradient:	A/B (0:100) to A/B (80:20) in 5min	
Flow Rate:	1mL/min	
Temperature:	22°C	
Injection:	20 μ L	
Detection:	UV @ 214nm	
Sample:	1. p-Aminophenol 2. Benzoic acid 3. Acetaminophen 4. Pseudoephedrine 5. Butyl paraben 6. Chlorphenarnamine 7. Diphenhydramine 8. Dextromethorphan	




LUNA CN (CYANO)


Proven Reproducibility

For carboxyl, carbonyl and amine containing compounds, Luna CN offers a unique polar selectivity in reversed phase and normal phase mode. Luna CN provides sharp peaks and great reproducibility from run-to-run, column-to-column and batch-to-batch. Luna's smooth silica allows for a more uniform bonding with improved resistance to bonded phase hydrolysis to produce one of the most stable CN phases. The result:


- Improved peak shape
- One of the most stable CN columns under reversed phase or normal phase conditions.
- pH stable from 1.5 to 7.0




Batch-to-Batch Reproducibility




Column-to-Column Reproducibility



Chromatographic Comparisons of CN Columns*



Phthalate Esters




* The comparative data presented here may not be representative for all applications.

LUNA NH₂ (AMINO)

Developed for Ruggedness

Luna NH₂ retains hydrogen-bonding compounds under three separation modes: Reversed Phased, Normal Phase, and Ion Exchange. Luna NH₂ provides reproducible retention and selectivity with improved column lifetime. Amino columns can be problematic as the bonded phase easily hydrolyses off the silica, shortening retention time over the life of the column. The bonded phase stability of Luna NH₂ is illustrated by 1.5 to 11.0 pH stability and 100% aqueous mobile phase stability. The result:

- Increased bonded phase stability improves reproducibility and column lifetime
- Excellent retention of simple sugars, complex sugars, sugar alcohols and other hydrogen-bonding compounds under Reversed Phase, Normal Phase, Ion Exchange conditions
- pH stable from 1.5 to 11.0
- Stable in 100% aqueous mobile phases




HPLC

Luna

Simple Sugars


App ID 14285
Column: Luna 5μ NH₂
Dimensions: 250 x 4.6mm
Order No: 00G-4378-E0

Mobile Phase: Acetonitrile/Water (80:20)
Flow Rate: 3mL/min
Temperature: 40°C
Detection: RI
Sample:
1. Xylose 5. Maltose
2. Fructose 6. Melezitose
3. Glucose 7. Raffinose
4. Sucrose




Steroids

App ID 14289
Column: Luna 5μ NH₂
Dimensions: 250 x 4.6mm
Order No: 00G-4378-E0
Mobile Phase: Hexane/Ethanol (85:15)
Flow Rate: 2mL/min
Temperature: 22°C
Detection: UV @ 240nm
Sample:
1. 11-Ketoprogesterone 4. Prednisolone 21-Acetate
2. 11-Hydroxyprogesterone 5. Cortisone
3. Cortisone Acetate 6. Prednisolone




12-18mer poly-DT Oligonucleotide

App ID 14628
Column: Luna 5μ NH₂
Dimensions: 150 x 4.6mm
Order No: 00F-4378-E0
Mobile Phase: A) 20mM Sodium Phosphate, 10% Acetonitrile, pH 7.0
B) 20mM Sodium Phosphate, 10% Acetonitrile,
1M Sodium Chloride, pH 7.0
Gradient: A/B (75:25) to A/B (5:95) in 50 minutes
Flow Rate: 1mL/min
Temperature: 40°C
Detection: UV @ 260nm
Sample: 12-18mer poly-DT oligonucleotide standard



Stability in 100% Aqueous Mobile Phase

Column: Luna 5μ NH₂
Dimensions: 250 x 4.6mm
Order No: 00G-4378-E0
Mobile Phase: 20mM Potassium Phosphate Buffer pH 2.7
Flow Rate: 1.5mL/min
Detector: UV@254nm
Temperature: Ambient
Injection: 2.5μL
Conditions: Continuously flushed at 1.0mL/min using 100% 20mM Potassium Phosphate Buffer pH 2.7 between injections



LUNA SCX (STRONG CATION EXCHANGE)


Luna SCX provides excellent resolution and peak shape of basic, cationic compounds. However, most SCX columns show poor peak shape and bad resolution causing many chromatographers to ignore this important phase for small molecule method development, until now. Luna SCX provides excellent durability over a wide pH range of 2 to 7.0 - wide pH stability is critical as mobile phase pH is a key factor in ion exchange chromatography. The result:

- **Resolving power and sharp peak shape to separate complex cationic/basic and nitrogen containing compounds**
- **5 and 10 μ columns and bulk media for analytical through preparative separations**
- **Good reproducibility assured by 25 QC tests and reported on a Materials Validation Document (MVD) accompanying each column**


App ID 14555

Childrens Tylenol Cold Syrup

Column: Luna 5 μ SCX 100Å
Dimension: 150 x 4.6mm
Order No: 00F-4398-E0
Mobile Phase: 50mM KH₂PO₄, pH 2.5/Acetonitrile (35:65)
Injection Volume: 1 μ L
Flow Rate: 1.5mL/min
Detection: UV @ 210nm
Sample Prep: Dissolve 1 part Childrens Tylenol Cold in 10 parts Methanol
Sample: 1. Acetaminophen
 2. Pseudoephedrine
 3. Dextromethorphan
 4. Chlorpheniramine




SCX Method Development and pH: The standard operating pH range for Luna SCX is 2.0 to 7.0. Most SCX methods are typically run between pH 2.0 and 5.0 for optimal performance. This ensures that nitrogen-containing analytes, especially those with adjacent conjugated system are protonated. Running in highly acidic (pH < 2.0) or basic (pH > 7.0) mobile phases may cause this phase to undergo degradation, as is common for all silica-based SCX phases.



App ID 14765

Peptides


Column: Luna 5 μ SCX
Dimensions: 150 x 4.6mm
Order No: 00F-4398-E0
Mobile Phase: A) 20mM Potassium Phosphate, 25% Acetonitrile, pH 2.5
 B) 20mM Potassium Phosphate, 25% Acetonitrile, 400mM Potassium Chloride, pH 2.5
Gradient: A/B (95:5) to A/B (10:90) in 45 minutes
Flow Rate: 1mL/min
Temperature: 35°C
Detection: UV @ 215nm
Injection Volume: 2 μ L (5 μ g on column)
Sample: Peptide Mixture - Substance P
 1. Fragment 5-11 (+1)
 2. Fragment 4-11 (+1)
 3. Fragment 2-11 (+2)
 4. Fragment 1-9 (+3)
 5. (+3)
 6. (ARG-PHE-TRP-LEU) (+3)



App ID 14767


Tryptic Digest of Bovine Cytochrome C

Column: Luna 5 μ SCX
Dimensions: 150 x 4.6mm
Order No: 00F-4398-E0
Mobile Phase: A) 20mM Potassium Phosphate, 25% Acetonitrile, pH 2.5
 B) 20mM Potassium Phosphate, 25% Acetonitrile, 350mM Potassium Chloride, pH 2.5
Gradient: 100% A to 100% B in 50 minutes
Flow Rate: 1mL/min
Temperature: 35°C
Detection: UV @ 215nm
Injection Volume: 50 μ L (20 μ g on column)
Sample: Bovine Cytochrome c trypsin digest



Luna SCX (continued)

Chromatographic Comparisons of SCX Columns*



Excellent Reproducibility

It is tough to meet Luna quality standards. The Silica Reproducibility and Bonded Phase Reproducibility Tanaka plots below illustrate the reproducibility of 13 different quality control parameters for 9 batches of Luna material - nearly no variation is observed. Each batch represents hours of quality control testing for over 25 different parameters. We know chromatographers rely on Luna for rugged, reproducible methods worldwide; therefore, an extensive Materials Validation Document (MVD) accompanying each column confirms consistent column specifications and consistent performance.

Material Validation Document



Quality Proven


Extensive testing and quality control test data is shown on the Materials Validation Document (MVD) that accompanies each 3 and 5 μ Luna column. Over 25 individual chemical and physical tests are performed and documented on the MVD, including:

- **Particle Analysis**
 - Particle Size
 - Pore Diameter
 - Particle Size Distribution
 - Surface Area
- **Total Metal Content**
- **Bonded Phase Coverage**
 - Total Carbon
 - Surface Coverage
- **Diagnostic Chromatography Test**
 - Inertness
 - Metal Sensitivity
 - Hydrophobic Index
- **Performance Chromatography Test**
 - Basic Drug Standards
 - Longevity Tests (pH Stability)
 - Scanning Electron Microscopy (measures surface smoothness and particle shape)



Column-to-Column Reproducibility

The chromatograms show consistency of inertness and hydrophobicity for Luna 5 μ C18(2) columns from 9 different batches. Almost no variation is observed.




HPLC

Luna


Silica Reproducibility

Quality control test data designed to monitor subtle changes in Luna silica that may affect reproducibility.



Bonded Phase Reproducibility

Bonding reproducibility test data from 9 Luna 5 μ C18(2) batches. Since differences in bonded phase coverage can show selectivity changes, each Luna batch undergoes 6 tests for bonding.



ORDERING INFORMATION

3 μ and 5 μ Capillary Columns (mm)

Phases	50 x 0.30	150 x 0.30	250 x 0.30	50 x 0.50	150 x 0.50	250 x 0.50
3 μ C8(2)	00B-4248-AC	00F-4248-AC	00G-4248-AC	00B-4248-AF	00F-4248-AF	00G-4248-AF
3 μ C18(2)	00B-4251-AC	00F-4251-AC	00G-4251-AC	00B-4251-AF	00F-4251-AF	00G-4251-AF
5 μ C8(2)	00B-4249-AC	00F-4249-AC	00G-4249-AC	00B-4249-AF	00F-4249-AF	00G-4249-AF
5 μ C18(2)	00B-4252-AC	00F-4252-AC	00G-4252-AC	00B-4252-AF	00F-4252-AF	00G-4252-AF



See p. 265 for Fused Silica Capillary Adapter and Capillary Guard Columns.



HPLC

Luna

MercuryMS™ LC/MS Cartridges (mm)

Qty Pricing	10 x 2.0	10 x 4.0	20 x 2.0	20 x 4.0
3 μ				
LUNA C18(2)	00N-4251-B0-CE	00N-4251-D0-CE	00M-4251-B0-CE	00M-4251-D0-CE
LUNA C8(2)	00N-4248-B0-CE	00N-4248-D0-CE	00M-4248-B0-CE	00M-4248-D0-CE
5 μ				
LUNA C18(2)	00N-4252-B0-CE	00N-4252-D0-CE	00M-4252-B0-CE	00M-4252-D0-CE
LUNA C8(2)	00N-4249-B0-CE	00N-4249-D0-CE	00M-4249-B0-CE	00M-4249-D0-CE



MercuryMS Cartridges require holder, see p. 268

SecurityGuard™ Cartridges require universal holder Order No.: KJ0-4282

3 μ Microbore and Minibore Columns (mm)

Phases	30 x 1.0	50 x 1.0	150 x 1.0	30 x 2.0	50 x 2.0	100 x 2.0	150 x 2.0	4 x 2.0mm /10pk	4 x 3.0mm /10pk
Silica(2)	—	—	—	00A-4162-B0	00B-4162-B0	00D-4162-B0	00F-4162-B0	AJ0-4347	AJ0-4348
C8(2)	00A-4248-A0	00B-4248-A0	00F-4248-A0	00A-4248-B0	00B-4248-B0	00D-4248-B0	00F-4248-B0	AJ0-4289	AJ0-4290
C18(2)	00A-4251-A0	00B-4251-A0	00F-4251-A0	00A-4251-B0	00B-4251-B0	00D-4251-B0	00F-4251-B0	AJ0-4286	AJ0-4287
CN	—	—	—	00A-4254-B0	00B-4254-B0	00D-4254-B0	00F-4254-B0	AJ0-4304	AJ0-4305
Phenyl-Hexyl	00A-4256-A0	00B-4256-A0	—	00A-4256-B0	00B-4256-B0	00D-4256-B0	00F-4256-B0	AJ0-4350	AJ0-4351
NH ₂	—	—	00F-4377-A0	00A-4377-B0	00B-4377-B0	00D-4377-B0	00F-4377-B0	AJ0-4301	AJ0-4302

for ID: 2.0-3.0mm 3.2-8.0mm

3 μ Analytical Columns (mm)

Phases	30 x 3.0	50 x 3.0	150 x 3.0	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	4 x 2.0mm /10pk	4 x 3.0mm /10pk
Silica(2)	—	—	00F-4162-Y0	—	00B-4162-E0	—	00D-4162-E0	00F-4162-E0	AJ0-4347	AJ0-4348
C8(2)	00A-4248-Y0	00B-4248-Y0	00F-4248-Y0	00A-4248-E0	00B-4248-E0	00C-4248-E0	00D-4248-E0	00F-4248-E0	AJ0-4289	AJ0-4290
C18(2)	00A-4251-Y0	00B-4251-Y0	00F-4251-Y0	00A-4251-E0	00B-4251-E0	00C-4251-E0	00D-4251-E0	00F-4251-E0	AJ0-4286	AJ0-4287
CN	—	00B-4254-Y0	00F-4254-Y0	00A-4254-E0	00B-4254-E0	00C-4254-E0	00D-4254-E0	00F-4254-E0	AJ0-4304	AJ0-4305
Phenyl-Hexyl	—	00B-4256-Y0	00F-4256-Y0	00A-4256-E0	00B-4256-E0	00C-4256-E0	00D-4256-E0	00F-4256-E0	AJ0-4350	AJ0-4351
NH ₂	—	00B-4377-Y0	00F-4377-Y0	00A-4377-E0	00B-4377-E0	—	00D-4377-E0	00F-4377-E0	AJ0-4301	AJ0-4302

for ID: 2.0-3.0mm 3.2-8.0mm



If Luna does not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, send in your comparative data within 45 days and keep the Luna column for FREE.

ORDERING INFORMATION

SecurityGuard™ Cartridges require universal holder Order No.: KJ0-4282

Phases	SecurityGuard™ Cartridges							
	50 x 1.0	150 x 1.0	250 x 1.0	30 x 2.0	50 x 2.0	150 x 2.0	250 x 2.0	4 x 2.0mm /10pk
Silica(2)	—	—	—	—	00B-4274-B0	00F-4274-B0	00G-4274-B0	AJ0-4347 AJ0-4348
C5	—	—	—	00A-4043-B0	00B-4043-B0	00F-4043-B0	—	AJ0-4292 AJ0-4293
C8	—	—	—	00A-4040-B0	00B-4040-B0	00F-4040-B0	—	AJ0-4289 AJ0-4290
C8 (2)	00B-4249-A0	00F-4249-A0	—	00A-4249-B0	00B-4249-B0	00F-4249-B0	00G-4249-B0	AJ0-4289 AJ0-4290
C18	—	—	—	00A-4041-B0	00B-4041-B0	00F-4041-B0	00G-4041-B0	AJ0-4286 AJ0-4287
C18 (2)	00B-4252-A0	00F-4252-A0	00G-4252-A0	00A-4252-B0	00B-4252-B0	00F-4252-B0	00G-4252-B0	AJ0-4286 AJ0-4287
CN	—	—	—	00A-4255-B0	00B-4255-B0	00F-4255-B0	00G-4255-B0	AJ0-4304 AJ0-4305
Phenyl-Hexyl	00B-4257-A0	00F-4257-A0	—	00A-4257-B0	00B-4257-B0	00F-4257-B0	00G-4257-B0	AJ0-4350 AJ0-4351
NH ₂	—	00F-4378-A0	—	00A-4378-B0	00B-4378-B0	00F-4378-B0	00G-4378-B0	AJ0-4301 AJ0-4302 for ID: 2.0-3.0mm 3.2-8.0mm



HPLC

Luna

Phases	SecurityGuard™ Cartridges							
	30 x 3.0	50 x 3.0	150 x 3.0	250 x 3.0	30 x 4.6	50 x 4.6	75 x 4.6	4 x 2.0mm /10pk
Silica(2)	—	—	—	—	—	00B-4274-E0	—	AJ0-4347 AJ0-4348
C5	—	—	00F-4043-Y0	—	—	00B-4043-E0	—	AJ0-4292 AJ0-4293
C8	—	—	00F-4040-Y0	00G-4040-Y0	00A-4040-E0	00B-4040-E0	00C-4040-E0	AJ0-4289 AJ0-4290
C8(2)	—	00B-4249-Y0	00F-4249-Y0	00G-4249-Y0	00A-4249-E0	00B-4249-E0	00C-4249-E0	AJ0-4289 AJ0-4290
C18	—	—	00F-4041-Y0	00G-4041-Y0	00A-4041-E0	00B-4041-E0	00C-4041-E0	AJ0-4286 AJ0-4287
C18(2)	00A-4252-Y0	00B-4252-Y0	00F-4252-Y0	00G-4252-Y0	00A-4252-E0	00B-4252-E0	00C-4252-E0	AJ0-4286 AJ0-4287
CN	—	00B-4255-Y0	00F-4255-Y0	00G-4255-Y0	00A-4255-E0	00B-4255-E0	00C-4255-E0	AJ0-4304 AJ0-4305
Phenyl-Hexyl	00A-4257-Y0	00B-4257-Y0	00F-4257-Y0	00G-4257-Y0	00A-4257-E0	00B-4257-E0	00C-4257-E0	AJ0-4350 AJ0-4351
NH ₂	—	00B-4378-Y0	00F-4378-Y0	00G-4378-Y0	00A-4378-E0	00B-4378-E0	00C-4378-E0	AJ0-4301 AJ0-4302 for ID: 2.0-3.0mm 3.2-8.0mm
SCX	—	—	—	—	—	00B-4398-E0	—	AJ0-4307 AJ0-4308





ORDERING INFORMATION

SecurityGuard™ Cartridges require universal holder Order No.: KJ0-4282

Phases	5μ Analytical, Semi-Prep and Preparative Columns (mm)								SecurityGuard™ Cartridges		
	100 x 4.6	150 x 4.6	250 x 4.6	250 x 10	250 x 15	250 x 21.2	250 x 30	4 x 2.0mm /10pk	4 x 3.0mm /10pk	10 x 10mm /3pk	
Silica(2)	00D-4274-E0	00F-4274-E0	00G-4274-E0	00G-4274-N0	—	00G-4274-P0	00G-4274-U0	AJ0-4347	AJ0-4348	AJ0-7223	
C5	00D-4043-E0	00F-4043-E0	00G-4043-E0	00G-4043-N0	—	—	—	AJ0-4292	AJ0-4293	AJ0-7372	
C8	00D-4040-E0	00F-4040-E0	00G-4040-E0	00G-4040-N0	—	—	—	AJ0-4289	AJ0-4290	AJ0-7222	
C8(2)	00D-4249-E0	00F-4249-E0	00G-4249-E0	00G-4249-N0	00G-4249-AK	00G-4249-P0	00G-4249-U0	AJ0-4289	AJ0-4290	AJ0-7222	
C18	00D-4041-E0	00F-4041-E0	00G-4041-E0	00G-4041-N0	—	—	—	AJ0-4286	AJ0-4287	AJ0-7221	
C18(2)	00D-4252-E0	00F-4252-E0	00G-4252-E0	00G-4252-N0	00G-4252-AK	00G-4252-P0	00G-4252-U0	AJ0-4286	AJ0-4287	AJ0-7221	
CN	00D-4255-E0	00F-4255-E0	00G-4255-E0	00G-4255-N0	—	—	—	AJ0-4304	AJ0-4305	AJ0-7313	
Phenyl-Hexyl	00D-4257-E0	00F-4257-E0	00G-4257-E0	00G-4257-N0	00G-4257-AK	00G-4257-P0	00G-4257-U0	AJ0-4350	AJ0-4351	AJ0-7314	
NH ₂	00D-4378-E0	00F-4378-E0	00G-4378-E0	00G-4378-N0	—	00G-4378-P0	00G-4378-U0	AJ0-4301	AJ0-4302	AJ0-7364	
SCX	00D-4398-E0	00F-4398-E0	00G-4398-E0	00G-4398-N0	—	00G-4398-P0	—	AJ0-4307	AJ0-4308	AJ0-7369	

for ID: 2.0-3.0mm 3.2-8.0mm 9-16mm

Phases	10μ Semi-Prep and Preparative Columns (mm)								SecurityGuard™ Cartridges	
	250 x 4.6	250 x 10	250 x 15	250 x 21.2	250 x 30	250 x 50	50 x 50	10 x 10mm /3pk	Guards	Cartridges
Silica(2)	00G-4091-E0	00G-4091-N0	—	00G-4091-P0	—	00G-4091-V0	03B-4091-V0	AJ0-7223		
C5	00G-4092-E0	00G-4092-N0	—	00G-4092-P0	—	—	—	AJ0-7372		
C8	00G-4093-E0	00G-4093-N0	—	00G-4093-P0	—	—	—	AJ0-7222		
C8(2)	00G-4250-E0	00G-4250-N0	00G-4250-AK	00G-4250-P0	00G-4250-U0	00G-4250-V0	03B-4250-V0	AJ0-7222		
C18	00G-4094-E0	00G-4094-N0	—	00G-4094-P0	—	—	—	AJ0-7221		
C18(2)	00G-4253-E0	00G-4253-N0	00G-4253-AK	00G-4253-P0	00G-4253-U0	00G-4253-V0	03B-4253-V0	AJ0-7221		
CN	00G-4300-E0	00G-4300-N0	—	00G-4300-P0	—	00G-4300-V0	03B-4300-V0	AJ0-7313		
Phenyl-Hexyl	00G-4285-E0	00G-4285-N0	—	00G-4285-P0	00G-4285-U0	00G-4285-V0	03B-4285-V0	AJ0-7314		
NH ₂	00G-4379-E0	00G-4379-N0	—	00G-4379-P0	—	—	—	AJ0-7364		
SCX	00G-4401-E0	00G-4401-N0	—	00G-4401-P0	—	—	—	AJ0-7369		

for ID: 9-16mm

Phases	15μ Preparative and Pilot Scale Columns (mm)					Guards	
	250 x 4.6	250 x 21.2	250 x 30	250 x 50	50 x 50	Guards	
Silica(2)	00G-4271-E0	00G-4271-P0	00G-4271-U0	00G-4271-V0	03B-4271-V0		
C8(2)	00G-4272-E0	00G-4272-P0	00G-4272-U0	00G-4272-V0	03B-4272-V0		
C18(2)	00G-4273-E0	00G-4273-P0	00G-4273-U0	00G-4273-V0	03B-4273-V0		
Phenyl-Hexyl	00G-4286-E0	00G-4286-P0	00G-4286-U0	00G-4286-V0	03B-4286-V0		



See p. 213 for Semi-Prep SecurityGuard Cartridge Holders



Improve analyte sensitivity and reduce baseline noise with Strata SPE tubes and well plates.
See pp. 9-26 for more information



COMBI-HTS™ PURIFICATION COLUMNS

New, 21.2mm and 30mm ID columns provide exceptional loadability and durability for high-throughput purifications. Luna Combi-HTS columns utilize the same Luna material, but undergo additional QC at higher flow rates to ensure excellent performance under high pressure and high flow conditions.

ORDERING INFORMATION

5µ Luna Combi-HTS Purification Columns (mm)					
Phases	50 x 21.2	100 x 21.2	50 x 30	75 x 30	100 x 30
Combi-HTS C8(2)	00B-4384-P0	00D-4384-P0	00B-4384-U0	00C-4384-U0	00D-4384-U0
Combi-HTS C18(2)	00B-4383-P0	00D-4383-P0	00B-4383-U0	00C-4383-U0	00D-4383-U0

METHOD DEVELOPMENT KITS

Luna Method Development Kits combine the excellent performance of Luna HPLC columns with either 3 or 5 different selectivities. The column kits include a Luna Phenyl-Hexyl (moderate polarity) and your choice of a Luna CN (polar) or a C8(2) or C18(2) (non-polars).

ORDERING INFORMATION

Method Development Kits					
Order No.	Description	Dimensions (mm)	Unit	Price	
KH0-4760	Luna Method Development Kit (5 columns) 1 each of 5µ C18(2), C8(2), C5, CN, Phenyl-Hexyl phases	150 x 4.6	5/pk		
KH0-4354	Luna Rapid-MD Kit (3 columns) 1 each of 3µ C18(2), C8(2) and Phenyl-Hexyl phases	50 x 4.6	3/pk		
KH0-4761	Luna Selectivity Kit 1 (3 columns) 1 each of 5µ C18(2), Phenyl-Hexyl, CN phases	150 x 4.6	3/pk		
KH0-4762	Luna Selectivity Kit 2 (3 columns) 1 each of 5µ C8(2), Phenyl-Hexyl, CN phases	150 x 4.6	3/pk		

Method Validation Kits					
Order No.	Description	Dimensions (mm)*	Unit	Price	
KH0-4355	Luna 5µ C18 Method Validation Kit	150 x 4.6	3/pk		
KH0-4356	Luna 5µ C18(2) Method Validation Kit	150 x 4.6	3/pk		
KH0-4357	Luna 5µ C8 Method Validation Kit	150 x 4.6	3/pk		
KH0-4358	Luna 5µ C8(2) Method Validation Kit	150 x 4.6	3/pk		
KH0-4359	Luna 5µ Phenyl-Hexyl Method Validation Kit	150 x 4.6	3/pk		

*Other dimensions available on request.



See p. 270 for LC/MS Computer-Based Training Program.



If Luna does not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, send in your comparative data within 45 days and keep the Luna column for FREE.



HPLC

Luna

METHOD VALIDATION KITS

Three columns of the same phase from three different batches provide an excellent benchmark for method validation. Method reproducibility is easily and quickly assessed, clearly demonstrating chromatographic performance validation potential and long term compliance with Luna Method Validation Kits.



BULK MATERIAL

ORDERING INFORMATION

10 μ Bulk Packings

Phases	100g	1kg	5kg	10kg	50kg	100kg
Silica(2)	04G-4091	04K-4091	04L-4091	04M-4091	04N-4091	04P-4091
Bonded Silica						
C5	04G-4092	04K-4092	04L-4092	04M-4092	04N-4092	04P-4092
C8	04G-4093	04K-4093	04L-4093	04M-4093	04N-4093	04P-4093
C8(2)	04G-4250	04K-4250	04L-4250	04M-4250	04N-4250	04P-4250
C18	04G-4094	04K-4094	04L-4094	04M-4094	04N-4094	04P-4094
C18(2)	04G-4253	04K-4253	04L-4253	04M-4253	04N-4253	04P-4253
CN	04G-4300	04K-4300	04L-4300	04M-4300	04N-4300	04P-4300
NH ₂	04G-4379	04K-4379	04L-4379	04M-4379	04N-4379	04P-4379
SCX	04G-4401	04K-4401	04L-4401	04M-4401	04N-4401	04P-4401
Phenyl-Hexyl						
Phenyl-Hexyl	04G-4285	04K-4285	04L-4285	04M-4285	04N-4285	04P-4285

10 μ -PREP Bulk PackingsLower cost excellent performance alternative to LUNA 10 μ .

Phases	100g	1kg	5kg	10kg	50kg	100kg
Silica(2)	04G-4322	04K-4322	04L-4322	04M-4322	04N-4322	04P-4322
Bonded Silica						
C8(2)	04G-4323	04K-4323	04L-4323	04M-4323	04N-4323	04P-4323
C18(2)	04G-4324	04K-4324	04L-4324	04M-4324	04N-4324	04P-4324
Phenyl-Hexyl						
Phenyl-Hexyl	04G-4325	04K-4325	04L-4325	04M-4325	04N-4325	04P-4325

15 μ Bulk Packings

Phases	100g	1kg	5kg	10kg	50kg	100kg
Silica(2)	04G-4271	04K-4271	04L-4271	04M-4271	04N-4271	04P-4271
Bonded Silica						
C8(2)	04G-4272	04K-4272	04L-4272	04M-4272	04N-4272	04P-4272
C18(2)	04G-4273	04K-4273	04L-4273	04M-4273	04N-4273	04P-4273
Phenyl-Hexyl						
Phenyl-Hexyl	04G-4286	04K-4286	04L-4286	04M-4286	04N-4286	04P-4286

