




# HPLC Chromatogram Index

## COSMOSIL HILIC



## Selection guide of mobile phase

COMOSIL HILIC column generates retention and separation by hydrophilic interaction (mainly hydrogen bond) and anion-exchange. Refer to following recommendations to select an appropriate mobile phase condition.

### (1) The effect of organic solvent type and content

- In general, acetonitrile/water is used as mobile phase.
- Retention increases as water content in the mobile phase decreased. (Fig.1)
- Use acetonitrile content in the mobile phase within the range of 0-95% (in general 50-95%).
- Methanol/water generates shorter retention than acetonitrile/water. (Fig.2)
- Use only HPLC grade solvents.




Fig.1 The effect of acetonitrile content on retention

Column; COSMOSIL HILIC

Mobile phase; Acetonitrile/ 10mmol/l CH<sub>3</sub>COONH<sub>4</sub>




Fig.2 Difference of acetonitrile and methanol on retention

Column; COSMOSIL HILIC

Mobile phase; Organic solvent/ 10mmol/l CH<sub>3</sub>COONH<sub>4</sub> = 90/10

### (2) The effect of buffer pH

- Keep pH of the mobile phase within the range of 2-7.5.
- The buffer around neutrality generates larger retention. (Fig.3)




Fig.3 The effect of buffer pH on retention

Column; COSMOSIL HILIC

Mobile phase; Acetonitrile / 10mmol/l buffer = 90/10

### (3) The effect of salt type and concentration

- When analyze ionic compounds, add salts or buffers in the mobile phase.
- When mobile phase has high acetonitrile content, note dissolubility of the salt. The dissolubility of phosphate buffers used widely in reversed phase chromatography is low in acetonitrile. Therefore use of phosphate buffers is not recommended. Keep the concentration of acetonitrile under 70% if use a phosphate buffer. Check that the salt does not precipitate when mixed with acetonitril before use.
- Ammonium acetate or formic acid ammonium buffers are recommended because they are soluble in high t acetonitrile content.

- Use the buffer concentration within the range of 5 - 100mmol/l. Moreover, Check that the salt does not precipitate after mixing buffer and acetonitrile.
- High salt concentration inhibits ion exchange and decreases retention. (Fig.4)
- Run mobile phase through membrane filter (less than 0.45µm in pore size) before use.




Fig.4 The effect of salt concentration on retention

Column: COSMOSIL HILIC

Mobile phase; Acetonitrile / 10mmol/l CH<sub>3</sub>COONH<sub>4</sub> = 50/50

#### (4) Selection of mobile phase


Following are the recommended mobile phases for different compound types.

Neutral compounds	→ Acetonitrile / Water = 90/10
Basic compounds	→ Acetonitrile / 10mmol/l CH <sub>3</sub> COONH <sub>4</sub> = 90/10
Amphoteric compounds	→ Acetonitrile / 10mmol/l CH <sub>3</sub> COONH <sub>4</sub> = 70/30
Acidic compounds	→ Acetonitrile / 10mmol/l CH <sub>3</sub> COONH <sub>4</sub> = 50/50
	↓ not eluted
	Acetonitrile / 10mmol/l Phosphate buffer (pH7.0)= 50/50

#### (5) Improvement of peak shape

Try following if peak shape is tailing. The peak shape might improve.

- Add 5mmol/l EDTA to mobile phase.
- Change to citrate buffer. (i. e. 10 mmol/l citrate buffer pH7.0)



#### (6) Others


- Use scrupulously degassed mobile phase. Air bubbles generate detection noise and accelerate column deterioration.
- We recommend keeping the chromatography conditions constant, since frequent changes of mobile phase shorten column life.

# Index

Sample name	Cas No.	Page	Sample name	Cas No.	Page
<b>A</b> Acetamide	60-35-5	1	L-Homoserine	672-15-1	14
Acetazolamide	59-66-5	1	Hydantoic Acid	462-60-2	14
Acetrizoic Acid	85-36-9	1	Hydantoin	461-72-3	14
Acrylic Acid	79-10-7	1	Hydroxylamine-O-sulfonic Acid	2950-43-8	14
L- $\alpha$ -Alanine	56-41-7	1	cis-4-Hydroxy-D-proline	2584-71-6	14
$\beta$ -Alanine	107-95-9	1	L-Hydroxyproline	51-35-4	15
Allantoic Acid	99-16-1	2	N-Hydroxysuccinimide	6066-82-6	15
Allantoin	97-59-6	2	<b>I</b> Indigo carmine	860-22-0	15
p-Aminobenzamide	3858-83-1	2	D-Isoascorbic Acid	89-65-6	15
p-Aminobenzoic Acid	150-13-0	2	Isocinchomeric Acid	100-26-5	15
4-Amino-n-butyric Acid [GABA]	56-12-2	2	L-Isoleucine	73-32-5	15
6-Aminohexanoic Acid	60-32-2	2	Isonicotinic Acid	55-22-1	16
5-Aminolevulinic Acid	106-60-5	3	Isonicotinohydrazide	54-85-3	16
2-Aminopyridine	504-29-0	3	Isopropyl $\beta$ -D-1-thiogalactopyranoside (IPTG)	367-93-1	16
3-Aminopyridine	462-08-8	3	<b>K</b> Kojic Acid	501-30-4	16
5-Amino-1H-tetrazole	4418-61-5	3	<b>L</b> L-Leucine	61-90-5	16
3-Amino-1H-1,2,4-triazole	61-82-5	3	N-D-Leucyl-L-tyrosine	3303-29-5	16
5-Aminouracil	932-52-5	3	L-Lysine	56-87-1	17
Amphotericin B	1397-89-3	4	<b>M</b> Maleic Acid	110-16-7	17
Angiotensin I(Human)	484-42-4	4	L-(-)-Malic Acid	97-67-6	17
Angiotensin II(Human)	4474-91-3	4	Malonic Acid	141-82-2	17
L-Arginine	74-79-3	4	Mecobalamin	13422-55-4	17
L(+)-Ascorbic Acid	50-81-7	4	Metanilic Acid	121-47-1	17
L-Asparagine	70-47-3	4	L-Methionine	63-68-3	18
L-Aspartic Acid	56-84-8	5	6-Methyl-2-thiouracil	56-04-2	18
6-Azauracil	461-89-2	5	N-Methylglucamine	6284-40-8	18
Aztreonam	78110-38-0	5	N-Methylhydroxylamine	593-77-1	18
<b>B</b> Benzamide	618-39-3	5	Mucic Acid	526-99-8	18
Benzenesulfonic Acid	98-11-3	5	Murexide	3051-09-0	18
Benzoic Acid	65-85-0	5	<b>N</b> Nicotinamide	98-92-0	19
Bromoacetic Acid	79-08-3	6	Nicotinic Acid	59-67-6	19
<b>C</b> Cacotheline	561-20-6	6	L-Noradrenaline	51-41-2	19
Camostat	59721-28-7	6	DL-Norleucine	616-06-8	19
L-Carnitine	541-15-1	6	DL-Norvaline	760-78-1	19
Ceftriaxone	73384-59-5	6	<b>O</b> L-Ornithine	70-26-8	19
Chloroacetic Acid	79-11-8	6	Orotic Acid	65-86-1	20
Citrazinic Acid	99-11-6	7	Oxalic Acid	144-62-7	20
Creatine	57-00-1	7	Oxamic Acid	471-47-6	20
Creatinine	60-27-5	7	Oxytocin	50-56-6	20
Cyanoacetic Acid	372-09-8	7	<b>P</b> D-Pantothenic Acid	79-83-4	20
Cyanuric Acid	108-80-5	7	L-(-)-Phenylalanine	63-91-2	20
L-Cysteine	52-90-4	7	p-Phenylenediamine	106-50-3	21
L-(-)-Cystine	56-89-3	8	L-(+)- $\alpha$ -Phenylglycine	2935-35-5	21
Cytidine	65-46-3	8	Phosphocreatine	67-07-2	21
Cytosine	71-30-7	8	O-Phospho-L-serine	407-41-0	21
<b>D</b> 3,4-Diaminobenzoic Acid	619-05-6	8	Picolinic acid	98-98-6	21
3,5-Diaminobenzoic Acid	535-87-5	8	Pivalic Acid	75-98-9	21
2,4-Diaminophenol	95-86-3	8	Procaterol	72332-33-3	22
DL-2,6-Diaminopimelic Acid	583-93-7	9	L-Proline	147-85-3	22
DL-2,3-Diaminopropionic Acid	54897-59-5	9	Propionic Acid	79-09-4	22
Diatrizoic Acid	117-96-4	9	Pyruvic Acid	127-17-3	22
Dipicolinic acid	499-83-2	9	<b>R</b> Ribose-5-phosphate	4300-28-1	22
Dithiouracil	2001-93-6	9	<b>S</b> D-Saccharic Acid	87-73-0	22
L-DOPA	59-92-7	9	Sarcosine	107-97-1	23
Dopamine	51-61-6	10	Sebacic Acid	111-20-6	23
meso-Erythritol	149-32-6	10	L-Serine	56-45-1	23
<b>F</b> Famotidin	76824-35-6	10	Sinigrin	3952-98-5	23
Folic Acid	59-30-3	10	Succinic Acid	110-15-6	23
Folinic Acid	58-05-9	10	Sulbactam	68373-14-8	23
Formamide	75-12-7	10	Sulfanilic acid	121-57-3	24
D-Fructose-6-phosphate	643-13-0	11	<b>T</b> L-(+)-Tartaric Acid	87-69-4	24
Fuchsine, Acid [Rubin S]	3244-88-0	11	Taurine	107-35-7	24
Fumaric Acid	110-17-8	11	L-Theanine	3081-61-6	24
<b>G</b> Gluconic Acid	526-95-4	11	2-Thiobarbituric Acid	504-17-6	24
$\alpha$ -D-Glucose-1-phosphate	59-56-3	11	2-Thiouracil	141-90-2	24
D-Glucose-6-phosphate	56-73-5	11	L-Threonine	72-19-5	25
D-Glucuronic Acid	6556-12-3	12	Todralazine	14679-73-3	25
L-Glutamic Acid	56-86-0	12	Trichloroacetic Acid	76-03-9	25
L-Glutamine	56-85-9	12	Tris(hydroxymethyl)aminomethane	77-86-1	25
Glutaric Acid	110-94-1	12	L-Tryptophan	73-22-3	25
DL-Glyceric Acid	600-19-1	12	L-Tyrosine	60-18-4	25
Glycinamide	598-41-4	12	<b>U</b> Uracil	66-22-8	26
Glycine	56-40-6	13	Urea	57-13-6	26
Glycolic Acid	79-14-1	13	Uridine	58-96-8	26
Glycylglycine	556-50-3	13	<b>V</b> L-Valine	72-18-4	26
Guanidoacetic Acid	352-97-6	13			
<b>H</b> 1,2,6-Hexanetriol	106-69-4	13			
L-Histidine	71-00-1	13			
L-Homocystine	626-72-2	14			


### COSMOSIL Chromatogram Index

Sample: Acetamide  
 CAS No.: [60-35-5]  
 Molecular formula:  $C_2H_5NO$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/  $H_2O=95/5$   
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 4.75min  
 Capacity factor: 0.57




### COSMOSIL Chromatogram Index

Sample: Acetazolamide  
 CAS No.: [59-66-5]  
 Molecular formula:  $C_4H_5N_3O_2S_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.2mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.99min  
 Capacity factor: 1.05




### COSMOSIL Chromatogram Index

Sample: Acetrisoic Acid  
 CAS No.: [85-36-9]  
 Molecular formula:  $C_9H_7I_3NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.8mg/ml  
 Injection volume: 1.0µl  
 Retention time: 16.39min  
 Capacity factor: 4.76




### COSMOSIL Chromatogram Index

Sample: Acrylic Acid  
 CAS No.: [79-10-7]  
 Molecular formula:  $C_3H_4O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 15.05min  
 Capacity factor: 4.28




### COSMOSIL Chromatogram Index

Sample: L-α-Alanine  
 CAS No.: [56-41-7]  
 Molecular formula:  $C_3H_7NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 5.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 6.67min  
 Capacity factor: 153




### COSMOSIL Chromatogram Index

Sample: β-Alanine  
 CAS No.: [107-95-9]  
 Molecular formula:  $C_3H_7NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 7.38min  
 Capacity factor: 1.81




### COSMOSIL Chromatogram Index

Sample: Allantoic Acid  
 CAS No.: [99-16-1]  
 Molecular formula:  $C_2H_4N_4O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 5.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 4.45min  
 Capacity factor: 0.69




### COSMOSIL Chromatogram Index

Sample: Allantoin  
 CAS No.: [97-59-6]  
 Molecular formula:  $C_4H_6N_4O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 8.75min  
 Capacity factor: 2.02




### COSMOSIL Chromatogram Index

Sample: p-Aminobenzamidine  
 CAS No.: [3858-83-1]  
 Molecular formula:  $C_7H_9N_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.07min  
 Capacity factor: 1.31




### COSMOSIL Chromatogram Index

Sample: p-Aminobenzoic Acid  
 CAS No.: [150-13-0]  
 Molecular formula:  $C_7H_7NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.4mg/ml  
 Injection volume: 1.0µl  
 Retention time: 16.97min  
 Capacity factor: 4.91




### COSMOSIL Chromatogram Index

Sample: 4-Amino-n-butyric Acid  
 CAS No.: [56-12-2]  
 Molecular formula:  $C_4H_7NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 7.67min  
 Capacity factor: 1.92





### COSMOSIL Chromatogram Index

Sample: 6-Aminohexanoic Acid  
 CAS No.: [60-32-2]  
 Molecular formula:  $C_6H_{13}NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 7.98min  
 Capacity factor: 2.03



### COSMOSIL Chromatogram Index


Sample: 5-Aminolevulinic Acid  
 CAS No.: [5451-09-2]  
 Molecular formula:  $C_7H_9NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 12.80min  
 Capacity factor: 3.87



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index

Sample: 2-Aminopyridine  
 CAS No.: [504-29-0]  
 Molecular formula:  $C_5H_5N_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 1.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 4.25min  
 Capacity factor: 0.39



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 3-Aminopyridine  
 CAS No.: [462-08-8]  
 Molecular formula:  $C_5H_5N_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.1mg/ml  
 Injection volume: 1.0µl  
 Retention time: 4.05min  
 Capacity factor: 0.51



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 5-Amino-1H-tetrazole  
 CAS No.: [4418-61-5]  
 Molecular formula:  $CH_3N_5$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 15.76min  
 Capacity factor: 4.49



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 3-Amino-1H-1,2,4-triazole  
 CAS No.: [61-82-5]  
 Molecular formula:  $C_2H_3N_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.2mg/ml  
 Injection volume: 1.0µl  
 Retention time: 7.01min  
 Capacity factor: 1.42



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 5-Aminouracil  
 CAS No.: [932-52-5]  
 Molecular formula:  $C_4H_5N_3O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 7.01min  
 Capacity factor: 1.42



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Amphotericin B  
 CAS No.: [1397-89-3]  
 Molecular formula:  $C_{47}H_{73}NO_{17}$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.25mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.34min  
 Capacity factor: 0.99



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Angiotensin I(Human)  
 CAS No.: [484-42-4]  
 Molecular formula:  $C_{23}H_{32}N_4O_7$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.4mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.28min  
 Capacity factor: 0.84



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Angiotensin II(Human)  
 CAS No.: [4474-91-3]  
 Molecular formula:  $C_{50}H_{71}N_{13}O_{12}$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.4mg/ml  
 Injection volume: 0.5µl  
 Retention time: 6.56min  
 Capacity factor: 1.29



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: L- Arginine  
 CAS No.: [74-79-3]  
 Molecular formula:  $C_6H_{12}N_4O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 7.97min  
 Capacity factor: 1.95



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: L(+)-Ascorbic Acid  
 CAS No.: [50-81-7]  
 Molecular formula:  $C_6H_8O_6$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV245nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.2mg/ml  
 Injection volume: 3.0µl  
 Retention time: 17.80min  
 Capacity factor: 5.31



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index

Sample: L-Asparagine  
 CAS No.: [70-47-3]  
 Molecular formula:  $C_4H_8N_2O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 4.88min  
 Capacity factor: 0.80





NACALAI TESQUE, INC



### COSMOSIL Chromatogram Index


Sample: L-Aspartic Acid  
 CAS No.: [56-84-8]  
 Molecular formula:  $C_4H_7NO_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 19.79min  
 Capacity factor: 6.01



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: 6-Azauracil  
 CAS No.: [461-89-2]  
 Molecular formula:  $C_4H_4N_2O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.3mg/ml  
 Injection volume: 0.5µl  
 Retention time: 9.65min  
 Capacity factor: 2.19



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Aztreonam  
 CAS No.: [78110-38-0]  
 Molecular formula:  $C_{13}H_{17}N_5O_8S_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 2.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 17.57min  
 Capacity factor: 5.18



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Benzamidine  
 CAS No.: [618-39-3]  
 Molecular formula:  $C_7H_9N_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 7.16min  
 Capacity factor: 1.46



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Benzenesulfonic Acid  
 CAS No.: [98-11-3]  
 Molecular formula:  $C_6H_7O_2S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 11.54min  
 Capacity factor: 3.05



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Benzoic Acid  
 CAS No.: [65-85-0]  
 Molecular formula:  $C_7H_6O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 15.19min  
 Capacity factor: 4.29



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Bromoacetic Acid  
 CAS No.: [79-08-3]  
 Molecular formula:  $C_2H_3BrO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 15.04min  
 Capacity factor: 4.31



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Cacoetheline  
 CAS No.: [561-20-6]  
 Molecular formula:  $C_{21}H_{21}N_3O_7$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 12.19min  
 Capacity factor: 3.23



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Camostat  
 CAS No.: [59721-28-7]  
 Molecular formula:  $C_{20}H_{22}N_4O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV265 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 7.16min  
 Capacity factor: 1.47



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: L-Carnitine  
 CAS No.: [541-15-1]  
 Molecular formula:  $C_7H_{13}NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 2.0mg/ml  
 Injection volume: 1.5µl  
 Retention time: 6.96min  
 Capacity factor: 1.78



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Ceftriaxone  
 CAS No.: [73384-59-5]  
 Molecular formula:  $C_{18}H_{18}N_4O_5S_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 11.36min  
 Capacity factor: 3.05



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Chloroacetic Acid  
 CAS No.: [79-11-8]  
 Molecular formula:  $C_2H_2ClO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 14.69min  
 Capacity factor: 4.15



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Citrazinic Acid  
 CAS No.: [99-11-6]  
 Molecular formula:  $C_7H_5NO_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 14.16min  
 Capacity factor: 3.98



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Creatine  
 CAS No.: [57-00-1]  
 Molecular formula:  $C_4H_9N_3O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.35min  
 Capacity factor: 1.40



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Creatinine  
 CAS No.: [60-27-5]  
 Molecular formula:  $C_4H_7N_3O$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 8.93min  
 Capacity factor: 2.08



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Cyanoacetic Acid  
 CAS No.: [372-09-8]  
 Molecular formula:  $C_3H_3NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 11.78min  
 Capacity factor: 3.56



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Cyanuric Acid  
 CAS No.: [108-80-5]  
 Molecular formula:  $C_3H_3N_3O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 7.61min  
 Capacity factor: 1.68



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: L-Cysteine  
 CAS No.: [52-90-4]  
 Molecular formula:  $C_3H_7NO_2S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 2.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 4.69min  
 Capacity factor: 0.73



NACALAI TESQUE, INC


### COSMOSIL Chromatogram Index

Sample: L-(-)-Cystine  
 CAS No.: [56-89-3]  
 Molecular formula:  $C_4H_8N_2O_4S_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 6.42min  
 Capacity factor: 1.38




### COSMOSIL Chromatogram Index

Sample: Cytidine  
 CAS No.: [65-46-3]  
 Molecular formula:  $C_4H_7N_3O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 7.22min  
 Capacity factor: 1.58




### COSMOSIL Chromatogram Index

Sample: Cytosine  
 CAS No.: [71-30-7]  
 Molecular formula:  $C_4H_5N_3O$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 11.22min  
 Capacity factor: 2.87




### COSMOSIL Chromatogram Index

Sample: 3,4-Diaminobenzoic Acid  
 CAS No.: [619-05-6]  
 Molecular formula:  $C_7H_7N_2O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.10mg/ml  
 Injection volume: 4.0µl  
 Retention time: 16.13min  
 Capacity factor: 4.62




### COSMOSIL Chromatogram Index

Sample: 3,5-Diaminobenzoic Acid  
 CAS No.: [535-87-5]  
 Molecular formula:  $C_7H_7N_2O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.1mg/ml  
 Injection volume: 4.0µl  
 Retention time: 16.54min  
 Capacity factor: 4.76




### COSMOSIL Chromatogram Index

Sample: 2,4-Diaminophenol  
 CAS No.: [95-86-3]  
 Molecular formula:  $C_6H_7N_2O$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 4.40min  
 Capacity factor: 0.51




### COSMOSIL Chromatogram Index

Sample: DL-2,6-Diaminopimelic Acid  
 CAS No.: [583-93-7]  
 Molecular formula:  $C_7H_{14}N_2O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.5µl  
 Retention time: 6.93min  
 Capacity factor: 1.56




### COSMOSIL Chromatogram Index

Sample: DL-2,3-Diaminopropionic Acid  
 CAS No.: [54897-59-5]  
 Molecular formula:  $C_4H_8N_2O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 5.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 14.38min  
 Capacity factor: 4.52




### COSMOSIL Chromatogram Index

Sample: Diatrizoic Acid  
 CAS No.: [1117-96-4]  
 Molecular formula:  $C_{11}H_9I_3N_2O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.8mg/ml  
 Injection volume: 1.0µl  
 Retention time: 14.98min  
 Capacity factor: 4.26




### COSMOSIL Chromatogram Index

Sample: Dipicolinic acid  
 CAS No.: [499-83-2]  
 Molecular formula:  $C_7H_7NO_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.37min  
 Capacity factor: 1.23




### COSMOSIL Chromatogram Index

Sample: Dithiouracil  
 CAS No.: [2001-93-6]  
 Molecular formula:  $C_4H_4N_2S_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.2mg/ml  
 Injection volume: 1.5µl  
 Retention time: 15.60min  
 Capacity factor: 4.15




### COSMOSIL Chromatogram Index

Sample: L-DOPA  
 CAS No.: [59-92-7]  
 Molecular formula:  $C_9H_9NO_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 3.0mg/ml  
 Injection volume: 3.0µl  
 Retention time: 7.12min  
 Capacity factor: 1.72




### COSMOSIL Chromatogram Index

Sample: Dopamine  
 CAS No.: [51-61-6]  
 Molecular formula:  $C_8H_{11}NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 1.0mg/ml  
 Injection volume: 4.0µl  
 Retention time: 7.73min  
 Capacity factor: 1.96




### COSMOSIL Chromatogram Index

Sample: meso-Erythritol [meso-Erythrite]  
 CAS No.: [149-32-6]  
 Molecular formula:  $C_4H_{10}O_5$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/  $H_2O$ =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.78min  
 Capacity factor: 1.18




### COSMOSIL Chromatogram Index

Sample: Famotidin  
 CAS No.: [76824-35-6]  
 Molecular formula:  $C_{16}H_{15}N_7O_2S_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.25mg/ml  
 Injection volume: 2.0µl  
 Retention time: 5.99min  
 Capacity factor: 1.06




### COSMOSIL Chromatogram Index

Sample: Folic Acid  
 CAS No.: [59-30-3]  
 Molecular formula:  $C_{19}H_{19}N_7O_6$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.25mg/ml  
 Injection volume: 2.0µl  
 Retention time: 14.09min  
 Capacity factor: 3.95




### COSMOSIL Chromatogram Index

Sample: Folinic Acid  
 CAS No.: [58-05-9]  
 Molecular formula:  $C_{20}H_{23}N_7O_7$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.25mg/ml  
 Injection volume: 2.0µl  
 Retention time: 13.36min  
 Capacity factor: 3.68





### COSMOSIL Chromatogram Index

Sample: Formamide  
 CAS No.: [75-12-7]  
 Molecular formula:  $CH_3NO$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/  $H_2O$ =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 4.58min  
 Capacity factor: 0.52



### COSMOSIL Chromatogram Index



Sample: D-Fructose-6-phosphate  
 CAS No.: [643-13-0]  
 Molecular formula:  $C_6H_{12}O_9P$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI  
 Attenuation:  $4 \times 10^3$  RIU/FS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 5.0µl  
 Retention time: 12.16min  
 Capacity factor: 3.64



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Fuch sine, Acid  
 CAS No.: [3244-88-0]  
 Molecular formula:  $C_{20}H_{17}N_3Na_2O_9S_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.5µl  
 Retention time: 13.82min  
 Capacity factor: 3.85



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Fumaric Acid  
 CAS No.: [110-17-8]  
 Molecular formula:  $C_4H_4O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.2mg/ml  
 Injection volume: 0.5µl  
 Retention time: 10.63min  
 Capacity factor: 2.75



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Gluconic Acid  
 CAS No.: [526-95-4]  
 Molecular formula:  $C_6H_{12}O_7$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.15min  
 Capacity factor: 0.81



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: α-D-Glucose-1-phosphate  
 CAS No.: [59-56-3]  
 Molecular formula:  $C_6H_{13}O_9P$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI  
 Attenuation:  $4 \times 10^3$  RIU/FS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 5.0µl  
 Retention time: 12.26min  
 Capacity factor: 3.68



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: D-Glucose-6-phosphate  
 CAS No.: [56-73-5]  
 Molecular formula:  $C_6H_{13}O_9P$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI  
 Attenuation:  $4 \times 10^3$  RIU/FS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 5.0µl  
 Retention time: 12.95min  
 Capacity factor: 3.94



NACALAI TESQUE, INC



### COSMOSIL Chromatogram Index

Sample: D-Glucuronic Acid  
 CAS No.: [6556-12-3]  
 Molecular formula:  $C_6H_{10}O_7$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.45min  
 Capacity factor: 0.92




### COSMOSIL Chromatogram Index

Sample: L-Glutamic Acid  
 CAS No.: [56-86-0]  
 Molecular formula:  $C_5H_9NO_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 19.38min  
 Capacity factor: 5.87





### COSMOSIL Chromatogram Index

Sample: L-Glutamine  
 CAS No.: [56-85-9]  
 Molecular formula:  $C_5H_{10}N_2O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 7.50min  
 Capacity factor: 1.85




### COSMOSIL Chromatogram Index

Sample: Glutaric Acid  
 CAS No.: [110-94-1]  
 Molecular formula:  $C_5H_8O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 10.45min  
 Capacity factor: 2.68




### COSMOSIL Chromatogram Index

Sample: DL-Glyceric Acid  
 CAS No.: [600-19-1]  
 Molecular formula:  $C_3H_4O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 6.0mg/ml  
 Injection volume: 5.0µl  
 Retention time: 15.68min  
 Capacity factor: 4.50



### COSMOSIL Chromatogram Index


Sample: Glycinamide  
 CAS No.: [598-41-4]  
 Molecular formula:  $C_2H_4N_2O$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6,Atten=8  
 Sample conc.: 1.0mg/ml  
 Injection volume: 3.0µl  
 Retention time: 13.35min  
 Capacity factor: 3.64






### COSMOSIL Chromatogram Index

Sample: Glycine  
 CAS No.: [56-40-6]  
 Molecular formula:  $C_2H_5NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 5.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 7.29min  
 Capacity factor: 1.77




### COSMOSIL Chromatogram Index

Sample: Glycolic Acid  
 CAS No.: [79-14-1]  
 Molecular formula:  $C_2H_3O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 15.28min  
 Capacity factor: 4.39




### COSMOSIL Chromatogram Index

Sample: Glycylglycine  
 CAS No.: [556-50-3]  
 Molecular formula:  $C_4H_8N_2O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 1.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 6.40min  
 Capacity factor: 1.27




### COSMOSIL Chromatogram Index

Sample: Guanidoacetic Acid  
 CAS No.: [352-97-6]  
 Molecular formula:  $C_3H_5N_3O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.61min  
 Capacity factor: 1.51




### COSMOSIL Chromatogram Index

Sample: 1,2,6-Hexanetriol  
 CAS No.: [106-69-4]  
 Molecular formula:  $C_6H_{14}O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/  $H_2O$ =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 1.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 5.19min  
 Capacity factor: 0.80





### COSMOSIL Chromatogram Index

Sample: L-Histidine  
 CAS No.: [71-00-1]  
 Molecular formula:  $C_6H_9N_3O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.50mg/ml  
 Injection volume: 1.0µl  
 Retention time: 8.38min  
 Capacity factor: 2.19



### COSMOSIL Chromatogram Index



Sample: L-Homocystine  
 CAS No.: [626-72-2]  
 Molecular formula:  $C_4H_9N_2O_4S_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 2.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 13.41min  
 Capacity factor: 4.10



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: L-Homoserine  
 CAS No.: [672-15-1]  
 Molecular formula:  $C_4H_7NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 7.03min  
 Capacity factor: 1.67



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Hydantoinic Acid  
 CAS No.: [462-60-2]  
 Molecular formula:  $C_3H_5N_2O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 16.33min  
 Capacity factor: 4.72



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Hydantoin  
 CAS No.: [461-72-3]  
 Molecular formula:  $C_3H_4N_2O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.01min  
 Capacity factor: 0.66



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Hydroxylamine-O-sulfonic Acid  
 CAS No.: [2950-43-8]  
 Molecular formula:  $H_3NO_3S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 2.0mg/ml  
 Injection volume: 3.0µl  
 Retention time: 15.60min  
 Capacity factor: 5.24



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: cis-4-Hydroxy-D-proline  
 CAS No.: [2584-71-6]  
 Molecular formula:  $C_5H_7NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.96min  
 Capacity factor: 1.65



NACALAI TESQUE, INC


### COSMOSIL Chromatogram Index

Sample: L-Hydroxyproline  
 CAS No.: [51-35-4]  
 Molecular formula:  $C_5H_9NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.49min  
 Capacity factor: 1.47




### COSMOSIL Chromatogram Index

Sample: N-Hydroxysuccinimide  
 CAS No.: [6066-82-6]  
 Molecular formula:  $C_4H_5NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.1 mg/ml  
 Injection volume: 1.5µl  
 Retention time: 11.29min  
 Capacity factor: 3.22




### COSMOSIL Chromatogram Index

Sample: Indigo carmine  
 CAS No.: [860-22-0]  
 Molecular formula:  $C_{16}H_8N_2Na_2O_6S_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.2mg/ml  
 Injection volume: 1.0µl  
 Retention time: 7.82min  
 Capacity factor: 1.79




### COSMOSIL Chromatogram Index

Sample: D-Isoascorbic Acid  
 CAS No.: [89-65-6]  
 Molecular formula:  $C_6H_8O_6$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 245nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.2mg/ml  
 Injection volume: 3.0µl  
 Retention time: 17.26min  
 Capacity factor: 5.11




### COSMOSIL Chromatogram Index

Sample: Isocinchomeronic Acid  
 CAS No.: [100-26-5]  
 Molecular formula:  $C_7H_7NO_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.91min  
 Capacity factor: 1.07




### COSMOSIL Chromatogram Index

Sample: L-Isoleucine  
 CAS No.: [73-32-5]  
 Molecular formula:  $C_6H_{13}NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.12min  
 Capacity factor: 0.95



### COSMOSIL Chromatogram Index



Sample: Isonicotinic Acid  
 CAS No.: [55-22-1]  
 Molecular formula:  $C_6H_5NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 1.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 16.45min  
 Capacity factor: 4.78



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Isonicotinohydrazide  
 CAS No.: [54-85-3]  
 Molecular formula:  $C_6H_5N_3O$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV265 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.50mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.37min  
 Capacity factor: 0.85



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Isopropyl β-D-1-thiogalactopyranoside  
 CAS No.: [367-93-1]  
 Molecular formula:  $C_{19}H_{35}O_5S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O=90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 0.1mg/ml  
 Injection volume: 0.5µl  
 Retention time: 6.23min  
 Capacity factor: 1.15



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Kojic Acid  
 CAS No.: [501-30-4]  
 Molecular formula:  $C_6H_6O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV245 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.1mg/ml  
 Injection volume: 1.0µl  
 Retention time: 3.83min  
 Capacity factor: 0.46



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: L-Leucine  
 CAS No.: [61-90-5]  
 Molecular formula:  $C_6H_{13}NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 4.91min  
 Capacity factor: 0.87



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: D-Leucyl-L-tyrosine  
 CAS No.: [3303-29-5]  
 Molecular formula:  $C_{15}H_{22}N_2O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 254nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 7.79min  
 Capacity factor: 1.96



NACALAI TESQUE, INC


### COSMOSIL Chromatogram Index

Sample: L-Lysine  
 CAS No.: [56-87-1]  
 Molecular formula:  $C_6H_{14}N_2O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 8.46min  
 Capacity factor: 2.13




### COSMOSIL Chromatogram Index

Sample: Maleic Acid  
 CAS No.: [110-16-7]  
 Molecular formula:  $C_4H_4O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.1mg/ml  
 Injection volume: 0.5µl  
 Retention time: 7.45min  
 Capacity factor: 1.62




### COSMOSIL Chromatogram Index

Sample: L-(-)-Malic Acid  
 CAS No.: [97-67-6]  
 Molecular formula:  $C_4H_6O_5$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 10.55min  
 Capacity factor: 2.71




### COSMOSIL Chromatogram Index

Sample: Malonic Acid  
 CAS No.: [141-82-2]  
 Molecular formula:  $C_3H_4O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 10.78min  
 Capacity factor: 2.81




### COSMOSIL Chromatogram Index

Sample: Mecobalamin  
 CAS No.: [13422-55-4]  
 Molecular formula:  $C_{65}H_{91}CoN_{13}O_{14}P$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV266 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.22min  
 Capacity factor: 1.35




### COSMOSIL Chromatogram Index

Sample: Metanilic Acid  
 CAS No.: [121-47-1]  
 Molecular formula:  $C_6H_7NO_3S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 13.68min  
 Capacity factor: 3.80




### COSMOSIL Chromatogram Index

Sample: L-Methionine  
 CAS No.: [63-68-3]  
 Molecular formula:  $C_5H_{11}NO_2S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 4.15min  
 Capacity factor: 0.54




### COSMOSIL Chromatogram Index

Sample: 6-Methyl-2-thiouracil  
 CAS No.: [56-04-2]  
 Molecular formula:  $C_5H_6N_2OS$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.1mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.58min  
 Capacity factor: 0.84




### COSMOSIL Chromatogram Index

Sample: N-Methylglucamine  
 CAS No.: [6284-40-8]  
 Molecular formula:  $C_7H_{17}NO_5$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI  
 Attenuation:  $4 \times 10^{-5}$  RIU/FS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 15.52min  
 Capacity factor: 4.22




### COSMOSIL Chromatogram Index

Sample: N-Methylhydroxylamine  
 CAS No.: [593-77-1]  
 Molecular formula:  $CH_3NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 1.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 13.45min  
 Capacity factor: 4.21




### COSMOSIL Chromatogram Index

Sample: Mucic Acid  
 CAS No.: [526-99-8]  
 Molecular formula:  $C_6H_{10}O_8$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 10.27min  
 Capacity factor: 2.62




### COSMOSIL Chromatogram Index

Sample: Murexide  
 CAS No.: [3051-09-0]  
 Molecular formula:  $C_8H_8N_4O_6$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 13.47min  
 Capacity factor: 3.69




### COSMOSIL Chromatogram Index

Sample: Nicotinamide  
 CAS No.: [98-92-0]  
 Molecular formula:  $C_6H_6N_2O$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.2mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.40min  
 Capacity factor: 0.77




### COSMOSIL Chromatogram Index

Sample: Nicotinic Acid  
 CAS No.: [59-67-6]  
 Molecular formula:  $C_6H_5NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 16.67min  
 Capacity factor: 4.87




### COSMOSIL Chromatogram Index

Sample: L-Noradrenaline  
 CAS No.: [51-41-2]  
 Molecular formula:  $C_8H_{11}NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.47min  
 Capacity factor: 1.07




### COSMOSIL Chromatogram Index

Sample: DL-Norleucine  
 CAS No.: [616-06-8]  
 Molecular formula:  $C_6H_{13}NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 4.89min  
 Capacity factor: 0.86




### COSMOSIL Chromatogram Index

Sample: DL-Norvaline  
 CAS No.: [760-78-1]  
 Molecular formula:  $C_5H_{11}NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.43min  
 Capacity factor: 1.07





### COSMOSIL Chromatogram Index

Sample: L-Ornithine  
 CAS No.: [70-26-8]  
 Molecular formula:  $C_5H_{12}N_2O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 8.39min  
 Capacity factor: 2.10



### COSMOSIL Chromatogram Index


Sample: Orotic Acid  
 CAS No.: [65-86-1]  
 Molecular formula:  $C_4H_4N_2O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 1.0µl  
 Retention time: 15.24min  
 Capacity factor: 4.36



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Oxalic Acid  
 CAS No.: [144-62-7]  
 Molecular formula:  $C_2H_2O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 5.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 12.08min  
 Capacity factor: 3.27



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Oxamic Acid  
 CAS No.: [471-47-6]  
 Molecular formula:  $C_2H_3NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.1mg/ml  
 Injection volume: 1.0µl  
 Retention time: 4.83min  
 Capacity factor: 0.71



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Oxytocin  
 CAS No.: [50-56-6]  
 Molecular formula:  $C_{43}H_{66}N_{12}O_{12}S_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.4mg/ml  
 Injection volume: 0.5µl  
 Retention time: 3.71min  
 Capacity factor: 0.39



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: D-Pantothenic Acid  
 CAS No.: [79-83-4]  
 Molecular formula:  $C_9H_{17}NO_5$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 13.21min  
 Capacity factor: 3.60



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index

Sample: L-(-)-Phenylalanine  
 CAS No.: [63-91-2]  
 Molecular formula:  $C_9H_9NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 254nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.10min  
 Capacity factor: 0.94




NACALAI TESQUE, INC




### COSMOSIL Chromatogram Index

Sample: p-Phenylenediamine  
 CAS No.: [106-50-3]  
 Molecular formula: C<sub>6</sub>H<sub>8</sub>N<sub>2</sub>  
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.1mg/ml  
 Injection volume: 0.5µl  
 Retention time: 4.15min  
 Capacity factor: 0.36




### COSMOSIL Chromatogram Index

Sample: L-(+)-α-Phenylglycine  
 CAS No.: [2935-35-5]  
 Molecular formula: C<sub>8</sub>H<sub>9</sub>NO<sub>2</sub>  
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.96min  
 Capacity factor: 1.27




### COSMOSIL Chromatogram Index

Sample: Phosphocreatine  
 CAS No.: [67-07-2]  
 Molecular formula: C<sub>4</sub>H<sub>10</sub>N<sub>3</sub>O<sub>5</sub>P  
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 11.42min  
 Capacity factor: 3.00




### COSMOSIL Chromatogram Index

Sample: O-Phospho-L-serine  
 CAS No.: [407-41-0]  
 Molecular formula: C<sub>3</sub>H<sub>7</sub>NO<sub>6</sub>P  
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 3.0µl  
 Retention time: 9.19min  
 Capacity factor: 2.24




### COSMOSIL Chromatogram Index

Sample: Picolinic acid  
 CAS No.: [98-98-6]  
 Molecular formula: C<sub>6</sub>H<sub>7</sub>NO<sub>2</sub>  
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.03min  
 Capacity factor: 0.92





### COSMOSIL Chromatogram Index

Sample: Pivalic Acid  
 CAS No.: [75-98-9]  
 Molecular formula: C<sub>5</sub>H<sub>10</sub>O<sub>2</sub>  
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 12.14min  
 Capacity factor: 3.28




### COSMOSIL Chromatogram Index

Sample: Procaterol  
 CAS No.: [72332-33-3]  
 Molecular formula:  $C_{16}H_{21}N_2O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Ammonium acetate =85/15  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 6.17min  
 Capacity factor: 1.25




### COSMOSIL Chromatogram Index

Sample: L-Proline  
 CAS No.: [147-85-3]  
 Molecular formula:  $C_5H_9NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.83min  
 Capacity factor: 1.22




### COSMOSIL Chromatogram Index

Sample: Propionic Acid  
 CAS No.: [79-09-4]  
 Molecular formula:  $C_3H_6O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 14.85min  
 Capacity factor: 4.24





### COSMOSIL Chromatogram Index

Sample: Pyruvic Acid  
 CAS No.: [127-17-3]  
 Molecular formula:  $C_3H_4O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 1.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 3.97min  
 Capacity factor: 0.39





### COSMOSIL Chromatogram Index

Sample: Ribose-5-phosphate  
 CAS No.: [4300-28-1]  
 Molecular formula:  $C_5H_{11}O_8P$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI  
 Attenuation:  $4 \times 10^{-5}$  RIU/FS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 5.0µl  
 Retention time: 8.02min  
 Capacity factor: 2.06




### COSMOSIL Chromatogram Index

Sample: D-Saccharic Acid  
 CAS No.: [87-73-0]  
 Molecular formula:  $C_6H_{10}O_8$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 10.48min  
 Capacity factor: 2.69




### COSMOSIL Chromatogram Index

Sample: Sarcosine  
 CAS No.: [107-97-1]  
 Molecular formula:  $C_2H_5NO_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 210nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.30min  
 Capacity factor: 1.40




### COSMOSIL Chromatogram Index

Sample: Sebacic Acid  
 CAS No.: [111-20-6]  
 Molecular formula:  $C_{18}H_{34}O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.5µl  
 Retention time: 9.43min  
 Capacity factor: 2.28




### COSMOSIL Chromatogram Index

Sample: L-Serine  
 CAS No.: [56-45-1]  
 Molecular formula:  $C_3H_7NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 7.92min  
 Capacity factor: 2.01




### COSMOSIL Chromatogram Index

Sample: Sinigrin  
 CAS No.: [3952-98-5]  
 Molecular formula:  $C_{10}H_{16}KNO_9S_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 12.57min  
 Capacity factor: 3.38




### COSMOSIL Chromatogram Index

Sample: Succinic Acid  
 CAS No.: [110-15-6]  
 Molecular formula:  $C_4H_4O_4$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 10.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 10.64min  
 Capacity factor: 2.74





### COSMOSIL Chromatogram Index

Sample: Sulbactam  
 CAS No.: [68373-14-8]  
 Molecular formula:  $C_8H_{11}NO_5S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l  
 Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 aufs  
 Sample conc.: 5.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 10.86min  
 Capacity factor: 2.81



### COSMOSIL Chromatogram Index



Sample: Sulfanilic acid  
 CAS No.: [121-57-3]  
 Molecular formula:  $C_6H_7NO_3S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.10mg/ml  
 Injection volume: 1.0µl  
 Retention time: 13.87min  
 Capacity factor: 3.87



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: L-(+)-Tartaric Acid  
 CAS No.: [87-69-4]  
 Molecular formula:  $C_4H_6O_6$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 1.5µl  
 Retention time: 10.52min  
 Capacity factor: 2.70



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index



Sample: Taurine  
 CAS No.: [107-35-7]  
 Molecular formula:  $C_2H_7NO_3S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.40min  
 Capacity factor: 1.25



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: L-Theanine  
 CAS No.: [3081-61-6]  
 Molecular formula:  $C_7H_{14}N_2O_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 5.0mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.89min  
 Capacity factor: 1.21



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 2-Thiobarbituric Acid  
 CAS No.: [504-17-6]  
 Molecular formula:  $C_4H_4N_2O_2S$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.1mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.18min  
 Capacity factor: 0.82



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 2-Thiouracil  
 CAS No.: [141-90-2]  
 Molecular formula:  $C_4H_4N_2OS$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.1mg/ml  
 Injection volume: 0.5µl  
 Retention time: 6.38min  
 Capacity factor: 1.11



NACALAI TESQUE, INC


### COSMOSIL Chromatogram Index

Sample: L-Threonine  
 CAS No.: [72-19-5]  
 Molecular formula:  $C_4H_9NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 10.0mg/ml  
 Injection volume: 2.0µl  
 Retention time: 7.19min  
 Capacity factor: 1.73




### COSMOSIL Chromatogram Index

Sample: Todalrazine  
 CAS No.: [14679-73-3]  
 Molecular formula:  $C_{11}H_{12}N_4O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV240 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 4.56min  
 Capacity factor: 0.51




### COSMOSIL Chromatogram Index

Sample: Trichloroacetic Acid  
 CAS No.: [76-03-9]  
 Molecular formula:  $C_2HCl_3O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 1.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 11.83min  
 Capacity factor: 3.17




### COSMOSIL Chromatogram Index

Sample: Tris(hydroxymethyl)aminomethane  
 CAS No.: [77-86-1]  
 Molecular formula:  $C_4H_{11}NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD  
 Attenuation: Gain=6, Atten=8  
 Sample conc.: 2.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 6.47min  
 Capacity factor: 1.48




### COSMOSIL Chromatogram Index

Sample: L-Tryptophan  
 CAS No.: [73-22-3]  
 Molecular formula:  $C_{11}H_{12}N_2O_2$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0) =70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254 nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 0.5mg/ml  
 Injection volume: 0.5µl  
 Retention time: 5.14min  
 Capacity factor: 0.95




### COSMOSIL Chromatogram Index

Sample: L-Tyrosine  
 CAS No.: [60-18-4]  
 Molecular formula:  $C_9H_9NO_3$   
 Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV 254nm  
 Attenuation: 0.128 auFS  
 Sample conc.: 5.0mg/ml  
 Injection volume: 1.0µl  
 Retention time: 5.92min  
 Capacity factor: 1.25




### COSMOSIL Chromatogram Index

Sample: Uracil  
CAS No.: [66-22-8]  
Molecular formula:  $C_4H_4N_2O_2$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5µl  
Retention time: 5.33min  
Capacity factor: 0.84




### COSMOSIL Chromatogram Index

Sample: Urea  
CAS No.: [57-13-6]  
Molecular formula:  $CH_4N_2O$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/  $H_2O$ =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0µl  
Retention time: 6.12min  
Capacity factor: 1.15




### COSMOSIL Chromatogram Index

Sample: Uridine  
CAS No.: [58-96-8]  
Molecular formula:  $C_9H_{12}N_2O_6$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.1mg/ml  
Injection volume: 1.0µl  
Retention time: 8.30min  
Capacity factor: 1.86



### COSMOSIL Chromatogram Index

Sample: L-Valine  
CAS No.: [72-18-4]  
Molecular formula:  $C_6H_{11}NO_2$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 5.63min  
Capacity factor: 1.14



## ナカライテスク株式会社

〒604-0855 京都市中京区二条通烏丸西入東玉屋町498

Web site

<http://www.nacalai.co.jp>

試薬はここに

価格・納期のご照会

**0120-489-552**

製品に関する技術的なご照会

E-mail: [info-tech@nacalai.co.jp](mailto:info-tech@nacalai.co.jp)

TEL:075-211-2703 FAX:075-211-2673

■販売取扱店